


**GALIT**  
EXPERTS IN DENTISTRY

# DENTAL EQUIPMENT

25 years of Experience, Leadership, and Expert  
Solutions for Dentistry



# EXPERTS IN DENTISTRY

**GALIT** is a renowned medical equipment manufacturer based in Ukraine, known for its unique blend of high-quality European components, versatility, and ergonomic design.

## EXPERTISE ABOVE ALL

For over a quarter of a century, the core mission of the GALIT brand has remained unchanged: to develop reliable and durable dental equipment while providing innovative solutions for dental professionals worldwide. Our commitment to the durability of our products reflects our dedication to fostering long-lasting relationships with our customers.

## BRAND VALUES:

**Quality:** We adhere to international quality standards for medical devices.

**Trust:** Our actions are guided by respect, integrity, and sincerity.

**Care:** Our team and our customers inspire us to act with dedication and passion.

**Innovation:** We leverage the latest technologies to achieve optimal results.

**Partnership:** We collaborate with healthcare professionals and partners to succeed together.

**Dedication:** We work meticulously to deliver exceptional results and outstanding customer service.



GALIT OMNIPRATIQUE

## Aesthetics and Functionality

The **GALIT OMNIPRATIQUE** dental unit is a market leader in its segment, chosen by more than 4,000 top clinics worldwide for daily procedures. It combines ergonomics, advanced technology, and sophisticated design to meet the needs of modern dentistry.

### Key Advantages of GALIT OMNIPRATIQUE:

- The unit is constructed using only high-quality components from European manufacturers and features durable materials such as polyurethane, metal, and ceramics.
- It allows the selection of up to five instruments for the dentist console and up to four for the assistant console, enabling a full range of dental procedures.
- Customisation options include the type and colour of upholstery, patient chair programmes, the instrument delivery type on the dentist console, and a standard or touch control panel. These allow the unit to be configured to meet the doctor's specific needs and ensure maximum comfort for both the dentist and the patient.
- The unit's design provides service personnel and medical staff with easy access to components.





GALIT AUTONOME



## | Independence for a Comfortable Practice

The **GALIT AUTONOME** dental unit is the ideal solution for clinics where service lines cannot be installed, as it only requires a connection to the power supply.

### Key Advantages of GALIT AUTONOME:

- This completely independent dental unit does not need any connections to air, vacuum, water, or sewage networks.
- It features a hydraulic unit that includes a compressor with a dryer and a stainless-steel receiver, a suction device with a bactericidal filter, and both a water and a secretion tank.
- The stable design allows for the easy attachment of additional equipment, such as microscopes, X-ray units, or monitors for convenient work.
- Components from leading European manufacturers with warranty service.
- The unit operates at a low noise level, ensuring a comfortable atmosphere for the dentist and the patient.
- The bactericidal HEPA filter (H14) ensures the air in the dental office is protected from viruses and bacteria, filtering it with 99.9% efficiency.
- Customisation options include the model, patient chair colour, upholstery type, additional programmes, the instrument delivery type on the dentist console, and a standard or touch control panel.





GALIT LIGHT



## Minimalism That Keeps Up With the Times

The **GALIT LIGHT** dental unit features a minimalist design and an ergonomic, space-saving construction, making it a versatile solution for compact spaces.

### Key Advantages of GALIT LIGHT:

- Made with high-quality European components, ensuring durability and reliability.
- The brushless NUCLEUS LED micromotor by TKD features a Duopad control panel and operates at speeds ranging from 100 to 40,000 rpm. Designed for endodontic procedures, it offers Auto-Stop, Auto-Reverse, Auto-Forward, and Rocking/ Reciprocal functions.
- The Water-Air Supply Unit from Dürr Dental controls pressure and filters water and air, automatically closing the channels for enhanced safety.
- The three-function 3FP-S Mytronic is an ultra-lightweight syringe with a soft grip, ensuring safe handling and compatibility with disposable nozzles.
- The comfortable patient chair features electromechanical movements and customisable options.
- The integrated Mytronic clean water system ensures superior hygiene and safety.





GALIT PRO

## Mastery in the Details

The **GALIT PRO** dental unit is a premium model designed for dentists who demand the highest quality standards and appreciate sophisticated design. Reliable and functional, this PRO unit also adds a stylish touch to any modern clinic.

### Key Advantages of GALIT PRO:

- Its elegant Italian design highlights the clinic's premium status, creating a professional image that inspires patient trust and helps clinics stand out.
- It is constructed using only premium components from leading European manufacturers to ensure high performance, durability, and impeccable quality.
- Customisation options include the type and colour of the upholstery, patient chair programs, instrument delivery, and control panel type, allowing you to tailor the unit to your needs and ensure maximum patient comfort.
- The robust design allows for easy integration of additional equipment, such as microscopes, X-rays, or monitors, expanding capabilities and creating a fully equipped workspace.





**GALIT AMBIDEXTRE**



## Comfort, Aesthetics, and Functionality

Practical and versatile, the **GALIT AMBIDEXTRE** dental unit is the ideal solution for clinics serving both right- and left-handed dentists.

### Key Advantages of GALIT AMBIDEXTRE:

- A unique rotation mechanism allows the dentist and assistant consoles to be switched from the right to the left side in just a few seconds.
- You can choose from ceiling, floor, or chair mounting options for the light and dentist console.
- The five-instrument dentist console, with upper or lower instrument delivery, is equipped with five instruments and two control panels for convenient treatment process management.
- The assistant console features a swivel arm and independent channels for saliva and dust ejectors, allowing easy movement in all directions—right, left, up, and down. This design ensures comfortable operation for both the dentist and the assistant.
- Patient chairs are designed for maximum comfort, offering various upholstery options and additional programs.



GALIT CONSOLE

## | Ergonomics Shaped by Experience

The **GALIT CONSOLE** is a modern dental unit designed for dentists prioritising freedom of movement and comfort in their daily practice. It combines advanced technology and ergonomic design to ensure efficient dental work and patient comfort.

### Key Advantages of GALIT CONSOLE:

- The chair is attached to the hydraulic unit, providing more legroom for the dentist and assistant, allowing them to work closely with the patient.
- The dentist console features a modern design with an upper or lower hose delivery and five customisable instruments to suit specific needs.
- A standard or touchscreen control panel for chair and unit operation can be used, and the power of all dynamic instruments can be adjusted.
- The patient chair is designed for maximum comfort and offers additional features, such as back and seat adjustments, a zero-return program, foot control, and an emergency stop function.
- Components from leading European manufacturers with warranty service.
- A ceramic swivel spittoon with electronic controls for flushing and cup filling.





GALIT NO STRESS



## Discover New Possibilities for Organizing Your Workspace

**GALIT NO STRESS** is an ergonomic and practical solution for your dental practice. It allows both the dentist and assistant to easily adjust the location and configuration of the equipment, enabling quick adaptation to specific procedures and patient needs.

### Key Advantages of GALIT NO STRESS:

- The devices are completely independent of air, vacuum, and sewage connections.
- The concept of stand-alone devices is ideal for retrofitting existing clinics and very small offices.
- The equipment meets the highest safety and hygiene standards following European regulations. Modern bactericidal HEPA filters provide ultra-fine air purification, capturing microparticles as small as 0.3 microns.
- There is no need for utility pipes. The mobile, self-contained units can be easily moved between rooms within your dental practice.





GALIT CART-5



## The Expert Solution for Modern Dentists

The **GALIT CART-5** mobile dental unit is a compact, flexible solution for modern dentists. The mobile dentist console combines stylish design with versatility to create optimal conditions for effective patient treatment.

### Key Advantages of GALIT CART-5:

- It offers a choice of up to five instruments, tailored to different procedures.
- It comes with a standard or touch control panel, providing complete instrument control.
- It features one-handed height adjustment for the dentist console, utilising a pneumatic system for smooth movement.
- The compact design saves space and complements the modern style of the dental office.
- The wheeled base allows for easy repositioning of the unit during patient treatment..
- The integrated air preparation and power units maximize space efficiency.
- It includes a 2-liter clean water system for enhanced hygiene.
- The easy-to-remove holders for autoclavable instruments enhance hygiene.
- Quick and simple disinfection of the unit's surface.





GALIT CART AUTONOME



## The Optimal Combination of Autonomy and Flexibility

The **GALIT CART AUTONOME** portable dental unit operates without a utility connection, providing maximum convenience in daily dental practice.

### Key Advantages of GALIT CART AUTONOME:

- It features a built-in oil-free DK 50 compressor with a 3-litre stainless-steel receiver housed in a soundproof cabinet to minimise noise.
- It includes a 2-litre clean water system, enabling use without a direct water supply.
- The unit's mobility facilitates easy manoeuvring and use in any position relative to the patient.
- Locking casters ensure stability of the unit during procedures.
- Available with either a standard or touch control panel, providing comprehensive control over dynamic instruments.
- The simple and stylish design ensures easy cleaning and disinfection.
- Removable instrument holders, suitable for autoclaving, to ensure high standards of hygiene and sterility.





**GALIT CART AUTONOME WITH ASPIRATION**



## Work Without Boundaries or Compromises

The **GALIT CART AUTONOME WITH ASPIRATION** is a mobile dental unit that combines efficiency and convenience, even in tight spaces, thanks to its integrated autonomous systems

### Key Advantages of GALIT CART AUTONOME WITH ASPIRATION:

- It has a built-in oil-free DK 50 compressor with a 3-litre stainless-steel receiver to minimise noise.
- It has a suction device with a bactericidal filter to ensure a clean working environment.
- It features a 2-litre clean water system and a 5-litre secretion tank for autonomous operation.
- The compact design is ideal for use in confined spaces.
- The wheeled frame with locking casters ensures safety and ease of movement.
- It features three instruments for the dentist with customisable options and two ejectors (saliva and dust) for the assistant.
- The power control panel for all dynamic instruments ensures intuitive operation.
- The thoughtful design makes it easy to clean and disinfect the surface.





MOBI SUC EASY AND MOBI SUC COMBI




## High Mobility and Efficiency

**MOBI SUC EASY and MOBI SUC COMBI** — are designed to effectively remove fluid or air from the patient's mouth during treatment procedures.

### Key Advantages of MOBI SUC EASY and MOBI SUC COMBI:

- Thanks to its wheeled frame, the suction device can be easily moved throughout the dental office.
- Reduced noise levels during operation.
- Excellent suction power is achieved within just a few seconds.
- The maximum suction flow rate is 1100 L/min.
- Special anti-static casters with a built-in brake provide secure locking.
- The devices are completely autonomous – simply plug them in and start working.
- The bactericidal HEPA filter (H14) purifies the air in the dental office, removing 99.9% of viruses and bacteria.





## Innovative Solutions for Quality Education

**GALIT TABLES FOR TRAINING CENTERS** — modern equipment that provides students with hands-on dental experience and prepares them for future clinical practice.

- **GALIT Training Table with Monitor** – Equipped with essential tools, a training simulator, and a monitor for ergonomic and effective dental training.
- **GALIT Training Table with Microscope** – Features dental tools, a training simulator, and a built-in microscope for hands-on experience in real treatment conditions.
- **GALIT Training Table with Lamp** – Includes necessary tools, a training simulator, and a dental lamp to enhance practical skills for modern dental practice.

FREE YOUR SMILE



## Comfort, Aesthetics and Functionality

**GALIT Medical Furniture** — high-quality and functional furniture that ensures a comfortable working environment for medical professionals while meeting modern hygiene standards.

### Key Advantages of GALIT Medical Furniture:

- Made from durable 1 mm metal, ensuring structural stability and longevity.
- Each module is designed to optimize workspace efficiency and provide maximum convenience.
- Surfaces are easy to clean and comply with medical hygiene standards.
- Metal structures prevent the accumulation of bacteria and viruses.
- Furniture can be recycled at the end of its lifespan, minimizing environmental impact.
- A wide range of modules and 29 colors allows customized solutions for any medical facility.
- Stationary, mobile, floor-standing, and wall-mounted options ensure optimal space utilization.

# Impressive Quality. Inspiring Innovations



The GALIT team is involved in every project stage, from design to implementation. For custom installation options, the first step is the engineering design process, ensuring you receive a product supported by the manufacturer's guarantee.



Our production facility covers 5,000 square meters and is equipped with state-of-the-art machinery. We prioritise maximum accuracy, quality, and short lead times to deliver exceptional products.

With a professional engineering department, skilled production specialists, and partnerships with over 30 renowned dental equipment brands, GALIT crafts expert solutions tailored to meet the needs and preferences of even the most discerning customers.

We prioritise the durability of our dental units, as we value the long-term relationships we build with our customers.

## GALIT - EXPERTS IN DENTISTRY



# GALIT

## Manufacturer

PC Galit  
6E, 15 Kvitnya Str., Baikivtsi, Ternopil region,  
47711, UKRAINE

+380-800-502-998

office@galit.te.ua  
service@galit.te.ua

## European representative

Galdent CZ s.r.o.  
K Vapence 448, 41503 Teplice, Czech Republic

+420-775-471-378

export@galit.te.ua  
galdentcz@gmail.com



Version: 20250115

The technical specifications and appearance of the equipment are subject to change without prior notice